

7516429

ClearKlens Cleansinald SC Disinfectant - 25ml Unit Dose Format

Product Description

ClearKlens Cleansinald is a surface-active disinfectant based on a blend of quaternary ammonium and alkylamines.

The ClearKlens sterile range was developed to meet the requirements of the pharmaceutical industry.

Directions for Use

- ClearKlens Cleansinald SC is used at 0.5% v/v aqueous solution.
- Empty the contents of a 25ml bottle into 5L of suitable water.

Features

- 0.2µm filtered, aseptically filled and double bagged.
- Manufactured to GMP standards with full batch traceability.
- Supplied with Certificates of Analysis, Sterility and Irradiation.
- Broad spectrum of activity.
- Aseptically filled and double bagged.
- Highly compatible with materials commonly found in pharmaceutical cleanrooms.
- SDS available upon request.



Purchase Information

7516429

Packaging

50 x 25ml bottles per case

